

## Perfect ultrasonic vibration...

### NEWTRON BOOSTER



Non-LED Handpiece,  
Quick-water disconnect.  
(Delivered with a Newtron handpiece, a power supply,  
a N\*1 tip and an autoclavable universal wrench)

**BUY 1 x Newtron  
Booster (Ref: F60203),  
GET  
2 x Newtron Tips  
FREE**

### NEWTRON P5 B.LED



LED Handpiece, Quick-water disconnect.  
(Delivered with a Newtron Slim B.Led handpiece with blue ring, a  
white LED ring, a F.L.A.G kit (5x1.5ml), an Australian power supply  
and a Pure Newtron kit (n.1, n.1S, H3 tips and an autoclavable uni-  
versal wrench)

**BUY 1 x Newtron  
P5 B.LED (Ref: F61105),  
GET  
2 x Newtron Tips  
FREE**

### NEWTRON P5XS B.LED



LED Handpiece with integrated irrigation  
tank for any solution. (Delivered with a Newtron Slim  
B.Led handpiece with blue ring, a white LED ring, a F.L.A.G kit  
(5x1.5ml), an Australian power supply and a Pure Newtron kit (n.1,  
n.1S, H3 tips and an autoclavable universal wrench).

**BUY 1 x Newtron  
P5XS B.LED (Ref: F62105),  
GET  
2 x Newtron Tips  
FREE**

## NEWTRON® Tip Kits for all your procedures

- NEWTRON tips conserve teeth and meet all clinical requirements
- NEWTRON tips are developed by dentists for dentists. They will offer you specific solutions, comfort and visibility.

### EndoSuccess Apical Surgery Kit



### F00069

Contains five ultrasonic instruments  
AS3D - Ref: F00065, AS6D - Ref: F00079,  
AS9D - Ref: F00067, ASLD - Ref: F00080,  
ASRD - Ref: F00081

### Endo Success Canal Access Prep Kit



### F88180

Contains 1 autoclavable tip support  
+ 3 tips:  
(CAP1, CAP2 and CAP3 + 1 universal wrench)

### Endo Success Retreatment Kit



### F00737

Contains six ultrasonic instruments  
ET18D-F88017, ETBD-F88020, ET20 -F88011,  
ET25 - F88018, ET25S -F88021, ETPR -  
F88019, ET20D - F88013, ET25L - F88022  
SO4 -F88009

## Superior osseous surgery for **immediate implant placement**

### Less Trauma, Less Drama

With Piezotome® CUBE bone surgery can be performed much faster with maximum power output and patient safety. The CUBE is especially effective with atraumatic extractions. Socket preservation, decreased patient pain and post-operative swelling are enhanced with CUBE facilitating immediate implant placement

## PIEZOTOME

CUBE

**BUY 1 x PIEZOTOME® CUBE (F50108),  
GET 1 x Extraction II Tip Kit (F87546) FREE**



## IMPLANT CENTER 2

### The alliance of technologies for safe and atraumatic bone surgeries

Highly reliable I-SURGE LED motor with its high torque (6Ncm) and wide speed range (100-40000Rpm).

Boosted, the Piezotome® LED handpiece contains 6 ceramic rings for faster and secure surgeries.

Access to all dental treatments with the NEWTRON® LED handpiece.

**BUY 1 x Implant Center II® (F27168),  
GET 1 x Extraction II Tip Kit (F87546) FREE**



## Mini LED Standard – Simply Powerful

**Acteon Mini LED** Outperforms in power, efficiency and speed. The most powerful in its category with 1250mW/cm<sup>2</sup> (with optical guide diameter 7.5 mm), without temperature rise. With the optical guide diameter 5.5 mm the power can reach 2000mW/cm<sup>2</sup>.


**The Mini LED Fast.** Only 6 to 12 sec to polymerize 2mm of any composite

- **FAST MODE** The choice of rapidity. Full power during 3, 4, 5 or 10 sec, according to the composite thickness
- **PULSE MODE** Limits temperature rise. 5 or 10 flashes of 1 sec for a secure composite polymerization
- **RAMPING MODE** Avoids composite retraction. Progressive power increase during 9 or 15 sec to ensure perfect restoration sealing
- **SAFETY** Base-charger includes a radiometer for light power intensity check-up



**BUY 1 x Mini LED (Ref: F02621), GET 1 x 7.5mm Opalescent Light Guide FREE**



# NEWTRON TIPS



**Over 80 different tips, with exclusive and patented designs, alloys and coatings that won't damage the different surfaces treated**



**irrisafe**

(4 - Pack)



IRR 20-21mm - F43807  
IRR 20-25mm - F43808  
IRR 25-21mm - F43805  
IRR 25-25mm - F43806



ET40  
F88012



ET40D  
F88014

**endodontics**

(4 - Pack)

K-Files  
(4-Pack)



K10/21mm - F43710  
K10/25mm - F43712  
K15/21mm - F43715  
K15/25mm - F43717  
K25/21mm - F43725  
K25/25mm - F43727  
K30/21mm - F43730  
K20/25mm - F43732

**ENDO SUCCESS**

Retreatment



ET18D  
F88017



ETBD  
F88020



ET20  
F88011



ET25  
F88018



ET25S  
F88021



ET25L  
F88022



ETPR  
F88019

**ENDO SUCCESS**

Apical Surgery



AS3D  
F00065



AS6D  
F00079



AS9D  
F00067



ASLD  
F00080



ASRD  
F00081

**ENDO SUCCESS**

Canal Access Prep



CAP1  
F88181




CAP2  
F88182



CAP3  
F88183

**PerfectMARGIN**

Rounded



PM1  
F02250



PM2  
F02251



PM3  
F02252



PM4  
F02253

**PerfectMARGIN**

Shoulder



PMS1  
F02254



PMS2  
F02255



PMS3  
F02256

**PerfectMARGIN**


Veneers



PMV1  
F02021



PMV2  
F02022



PMV3  
F02023



PMV4  
F02024



PMV5  
F02025



PMV6  
F02026

ROUND TORQUE WRENCHES




F81322

**BUY 4 x Newtron tips**


**GET 1 x FREE** (mix N Match Ok)

## PHOTOGRAPHIC MIRRORS - PURE REFLECT® TECHNOLOGY

To meet modern dentistry standards, the exacting nature of the **ACTEON® PURE REFLECT®** technology offers true clarity with a reflection rate close to 100%, providing a precise image of colors and dental structures without distortion. Make the difference! Upgrade the quality of the dental pictures, enhance communication with patients and laboratories and highlight your clinical publications!



**BUY 3 x Photographic Mirrors**  
**GET 1 x Mirror Cassette (Ref # 233.19) FREE**




## MOUTH MIRRORS - PURE REFLECT® TECHNOLOGY

**ACTEON® PURE REFLECT®** Technology: Optimum solution for a perfect visual examination!



To meet modern dentistry standards, the exacting nature of the **ACTEON® PURE REFLECT®** mirrors offers a perfect image of the operating site with a reflection rate close to 100%, ensuring beyond the visible more reliable visual diagnosis.

Available on the complete range of mouth mirrors, photographic mirrors and micro-mirrors, **PURE REFLECT®** technology is your best allies from perfect diagnostic to accurate treatment, for more predictable results.

- 100% level of color reflexion
- 0 image distortion
- Image reflected directly on the surface of the glass, distortion free.
- No double image



For your **FREE Pure Reflect Mouth Mirror Sample** Size 5, call Presidential on **1800 773 743**, Email [Jamesw@presidential.com.au](mailto:Jamesw@presidential.com.au) OR [Luke@presidential.com.au](mailto:Luke@presidential.com.au)



**BUY 2 x Boxes of Acteon Mouth Mirrors** **GET 1 x Box FREE** (mix N Match Ok)

## Complete your order form and email to [orders@presidential.com.au](mailto:orders@presidential.com.au)

**BUY** 1 x Piezotome Cube (Ref# F87546)  
**GET** 1 x Extraction Tip kit (Ref# F87546) **FREE**

**BUY** 1 x Mini LED Standard  
**GET** 1 x 7.5mm Opalescent Light Guide **FREE**

**BUY** 1 x Implant 2 Centre (Ref# F27203)  
**GET** 1 x Extraction Tip kit (Ref# F87546) **FREE**

**BUY** 3 x Photographic Mirrors  
**GET** 1 x Cassette (Ref# 233.19) **FREE**

**BUY** 1 x Newtron Booster  
**GET** 2 x Newtron Tips **FREE**

**BUY** 2 x Boxes Mouth Mirrors  
**GET** 1 x Box **FREE** (mix N Match OK)

**BUY** 1 x Newtron P5 B.Led  
**GET** 2 x Newtron Tips **FREE**

For your **FREE Pure Reflect Mouth Mirror Sample**  
Size 5, call Presidential on **1800 773 743**,  
Email [Jamesw@presidential.com.au](mailto:Jamesw@presidential.com.au)  
OR [Luke@presidential.com.au](mailto:Luke@presidential.com.au)

**BUY** 1 x Newtron P5 XS Bled  
**GET** 2 x Newtron Tips **FREE**

**BUY** 4 x Newtron Tips  
**GET** 1 **FREE**

Book an Acteon Lunch'n Learn  
to experience the  
**Piezotome CUBE (LED)**  
with your Acteon Clinical Specialist and  
**Receive 1 Box Riskontrol & Adapter FREE**

For bookings call Luke or James on **1800 773 743**,  
Email [Jamesw@presidential.com.au](mailto:Jamesw@presidential.com.au)  
OR [Luke@presidential.com.au](mailto:Luke@presidential.com.au)

Please print **PROMO CODE: AUS N N2019-3** legibly on your Invoice. To process your free good(s), To process your free goods, email and attach a PDF of your proof of purchase to [acteon.anz.promo@acteongroup.com](mailto:acteon.anz.promo@acteongroup.com)

Name: \_\_\_\_\_

Office Name: \_\_\_\_\_

Address: \_\_\_\_\_


Postcode: \_\_\_\_\_

Phone: \_\_\_\_\_

Email: \_\_\_\_\_

Preferred Dealer: \_\_\_\_\_

Dealer Account #: \_\_\_\_\_



**presidential**  
dental | medical | veterinary

[www.presidential.com.au](http://www.presidential.com.au)